



EXCELLENCE IN ELECTROMAGNETIC THERAPY

Leading Rehabilitation Technology

ايفوميديكال
EVO MEDICAL
MEDICAL EQUIPMENT للمعدات الطبية

Dubai, UAE

Phone: +(971) 55 6401016

Website: www.evomedicals.com

Email: info@evomedicals.com





REMEDI

Since its foundation in 2003, REMEDI has tried to care for the patients suffering from different types of physical pains all over the world by developing world's 1st electromagnetic stimulation device.

Now experience the powerful values of electromagnetic stimulation devices from REMEDI.

- ✓ Highly convenient and comfortable treatment for users & patients: non-invasive & in-dress treatment, short time and intensive treatment session available
- ✓ Deep penetrating treatment effect by 3 Tesla electromagnetic field for more than 10 cm deep inside body: deep penetration into skin, muscle and bone as well
- ✓ Unique oil circulation cooling transducer for non-tiring and intensive treatment activities: technically proven life long transducer due to non-vaporization of insulating oil








REMEDI

All rights reserved.


All products and their specifications contained in the catalogue may differ from country to country where they are sold. Specification change may occur without prior notice.

DEEP CLINICAL EFFECTS

-  Fast pain relief either at chronic or acute stage
-  Joint blockage release
-  Enhancement of fracture healing process
-  Muscle relaxation
-  Strong effect for genitourinary disorders: urinary incontinence, pelvic pain


A patient with L4 / 5 HNP

BEFORE TREATMENT



Right leg & the sole of the right foot show low temperature due to pain and poor blood circulation.

AFTER TREATMENT



However, after some electromagnetic treatment, the treated areas show normal temperature.

GENERAL CLINICAL INDICATIONS



Optimus Pro

World's 1st Combo of
Electromagnetic & RSWT Solution



FEATURES & BENEFITS

- Unique and innovative combination of 3 Tesla electromagnetic stimulator & RSWT
- Perfect solution for acute and chronic musculoskeletal pain treatment
- User-friendly interface with 8" touch screen LCD display
- Pre-set operation mode : 60 kinds of operation mode
- Unique oil circulation cooling transducer for continuous treatment

VARIOUS APPLICATORS



※ Applicators are for RSWT only.

SPECIFICATION

Magnetic Power	3 Tesla
Magnetic Frequency	1~100 Hz
Magnetic Treatment Mode	Auto Mode : 20 Manual Mode : 20 User Mode : 20
Magnetic Treatment Time	1~60 min
RSWT Strength	1~6 Bar, 0.1 Bar Step
RSWT Frequency	1~21 Hz
RSWT Treatment mode	M1~M20
RSWT Treatment Shot	100~5000, continuous, 100 Step
RSWT Applicator	6, 10, 15, 20 & 35 mm (Option)
Interface	8" LCD Touch Screen, Jog Shuttle
Dimension	45 x 60 x 110 Cm
Weight	68.5 Kg

MAJOR CLINICAL APPLICATION



Shoulder



Back



Elbow

FIELD OF APPLICATION

- Orthopedics
- Rehabilitation
- Genitourologist
- Sports medicine



ايفوميديكال
EVO MEDICAL
MEDICAL EQUIPMENT للمعدات الطبية

Dubai, UAE

Phone: +(971) 55 6401016

Website: www.evomedicals.com

Email: info@evomedicals.com

