


Salus FSWT

Powerful, precise and well-defined FSWT System



REMED
Rehabilitation Medical Company

ايفوميديكال
EVO MEDICAL
MEDICAL EQUIPMENT للمعدات الطبية

Dubai, UAE

Website: www.evomedicals.com

Phone: + (971) 55 6401016

Email: info@evomedicals.com


FEATURES & BENEFITS

- Direct focusing technology by piezo-electric elements for minimized pain
- User-friendly interface on 8 inch touch screen
- More durable applicator
- Easy adjustment of penetration depth and energy levels
- Wide range of energy settings up to 12 Hz
- Plug & Play applicator

SPECIFICATION

Energy level	max 0.5 mJ/mm ² , 0~100%, 1% step
Frequency	1~12 Hz, 1 Hz step
Shocks	0~10,000 shocks, 100 shock step
Power	220V AC
Dimension	424X405X207 (mm)
Weight	38kg (With Trolley)

VARIOUS GEL PADS



(8 KINDS OF GEL PADS : 5, 10, 15, 20, 25, 30, 35, 40mm)

CLINICAL INDICATIONS

- Tennis elbow/ Golf elbow
- Calcific tendonitis : calcified shoulder
- Plantar talagia
- Degenerative arthritis
- Patella tis syndrome
- Greater trochanteric pain syndrome
- Erectile dysfunction

MAJOR CLINICAL APPLICATION

